

Hello Bio, Inc.  
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500  
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



---

## DATASHEET

### TBS-T with Tween 20 (20x)

---

#### Product overview

|                               |   |
|-------------------------------|---|
| <b>Name</b>                   | TBS-T with Tween 20 (20x)   |
| <b>Cat No</b>                 | HB6971  |
| <b>Biological description</b> | 20x Tris Buffered Saline wash buffer for use during immunolabeling techniques such as IHC. When diluted to 1x contains 130mM NaCl, 20mM Tris, 0.1% Tween 20, 0.002% ProClin-300, pH7.4. |
| <b>Biological action</b>      | Buffer  |
| <b>Description</b>            | 20x Tris Buffered Saline wash buffer to remove non-specific binding and for use during immunolabeling techniques such as IHC.   |

---

#### Biological Data

|                          |  |
|--------------------------|--|
| <b>Application notes</b> | This is a 20x buffer: to use, dilute in dH <sub>2</sub> O (50ml of 20x buffer + 950ml dH <sub>2</sub> O) |
|--------------------------|--|

---

#### Solubility & Handling

|                             |   |
|-----------------------------|---|
| <b>Storage instructions</b> | Room temperature  |
| <b>Important</b>            | This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use. |

---

#### Chemical Data

|                   |                    |
|-------------------|--------------------|
| <b>Appearance</b> | Colorless solution |
|-------------------|--------------------|

---