Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

Biotin Tyramide Signal Amplification Kit - Janelia Fluor® 525

Product overview

Name Cat No Biological description Biotin Tyramide Signal Amplification Kit - Janelia Fluor® 525 HB19791

Biotin Tyramide Signal Amplification Kit – Janelia Fluor® 525 enables immunofluorescent detection of low abundance targets with a bright and photostable fluorophore. Biotin Tyramide Signal Amplification Kits enhance signal through deposition of biotin tyramide on neighbouring molecules. Janelia Fluor® 525 detects biotin with extremely high affinity, resulting in increased sensitivity and signal detection. Janelia Fluor® 525 and the other members of the Janelia Fluor® family are bright and highly photostable fluorophores particularly suited for super resolution imaging such as dSTORM and STED.

This kit provides the reagents for 100 x 1mL reactions, with a protocol that can easily be adapted for immunohistochemistry (IHC), immunocytochemistry (ICC) and in situ hybridisation (ISH).

Key features:

- Janelia Fluor® 525 Streptavidin for detection (Ex: 534nm, Em: 559nm)
- Bright and photostable signal for repeated imaging
- Enhanced signal of hard-to-detect targets
- Reduced primary antibody use
- Suited for IHC(IF), ICC, FISH
- Suitable for super resolution imaging


fluorescence imaging, ICC, IF, IHC(IF)


Immunofluorescent signal amplification kit using biotin tyramide and Streptavidin Janelia Fluor® 525.


Images


Applications

Description









Biological Data

Application notes

Biotin Tyramide Signal Amplification Protocol

Please follow this link to the amplification protocol.

Solubility & Handling

Storage instructions Reconstitution advice

-20°C