

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

F-127 (20% solution in DMSO)

Product overview

Name	F-127 (20% solution in DMSO)
Cat No	HB9631
Biological description	Nonionic detergent commonly used to solubilize water-insoluble dyes. F 127 (20% solution in DMSO) is used to disperse AM (acetoxymethyl) esters of ion indicators including Fluo-4 AM UltraPure , Fura-2 AM (Cell permeant) and BAPTA-AM , as well as cell tracer dyes such as CFDA-SE. F-127 (10% solution in water, sterile-filtered) also available .
Alternative names	F-127; F127
Biological action	Reagent
Description	Nonionic detergent for solubilizing ion indicators (20% solution in DMSO)

Solubility & Handling

Storage instructions	Room temperature
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use

Chemical Data

CAS Number	9003-11-6
PubChem identifier	24751

References

Calcium rises abruptly and briefly throughout the cell at the onset of anaphase

Poenie et al (1986) Science.
