Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

(-)-Bicuculline methiodide

Product overview

Name (-)-Bicuculline methiodide

Cat No HB0893
Alternative names BIC, BMI
Biological action Antagonist
Purity >98%

Customer commentsGood quality product: (-)-Bicuculline methiodide is used routinely in our lab for a number of

experiments. It is shipped quickly, packaged well, dissolves without problem, and blocks GABAA-

receptor activity as it should! Verified customer, Sickkids foundation


We routinely use this compound from Hello Bio to inhibit GABA-A receptors in electrophysiological

recordings from rodent brain slices. Verified customer, University of Montana

Description Prototypic, competitive GABA_A receptor antagonist

Images







Biological Data

Biological description

Methiodide salt form of (+)-bicuculline.

Prototypic, competitive $GABA_A$ receptor antagonist which displaces GABA from the agonist binding site to prevent receptor activation.

Also acts as a negative allosteric inhibitor of channel opening to inhibit $GABA_A$ receptor activation by anaesthetic agents.

Additionally shows activity at SK calcium-activated potassium channels, nicotinic acetylcholine receptors and acetylcholinesterase.

Reversibly and competitively blocks GABA_A receptor mediated currents. Widely used to isolate glutamate receptor mediated EPSCs (excitatory postsynaptic potentials).

Shows convulsant action and induces epilepsy.

Freebase, methochloride and methobromide salts also available.

the actions of the neurotransmitter GABA. Bicuculline is commonly used at concentrations of 100 µM and above.

Bicuculline methiodide from Hello Bio reduces both spontaneous inhibitory post synaptic currents (IPSC) and evoked IPSCs (see Fig 1 above). It was effective at concentrations of 1 mM with complete receptor blockade at 100 µM.

#Protocol 1: Evoked and spontaneous inhibitory post synaptic currents (IPSCs)

- Whole cell voltage clamp recordings were obtained from layer V neurons of the mouse prelimbic cortex brain slice.
- A stimulating electrode was placed in layers II/III and IPSCs were evoked by a single square (150 µs) pulse every 10 sec at a stimulus intensity that gave a reliable IPSC.
- IPSCs were evoked at a range of neuron holding voltages to measure the reversal potential of the current to ensure it was GABAergic.
- Neurons were held at 0mV and IPSCs continuously stimulated and recorded in response to 5 min applications of varying concentrations of Bicuculline methiodide until complete receptor inhibition.
- Spontaneous IPSCs were recorded before and after addition of Bicuculline methiodide by holding the neuron at 0mV and recording for 10 sec.
- All recordings for IPSCs were made in the presence of AMPAR antagonists.

Solubility & Handling

Storage instructions Solubility overview

Handling

Important

Room temperature

Soluble in water (20mM) or DMSO (50mM)

This compound is light sensitive; exposure to light may affect compound performance. We therefore recommend storing the solid material and any solutions in the dark and protecting from light.


This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

Chemical Data

Chemical name [R-(R*.S*)]-5-(6.8-Dihvdro-8-oxofuro[3.4-e

Molecular Weight **Chemical structure** [-1,3-benzodioxol-6-yl)-5,6,7,8-tetrahydro-6,6-dimethyl-1,3-dioxolo[4,5-g]isoquinolinium iodide 509.3



Molecular Formula C21H20INO6 **CAS Number** 40709-69-1 **PubChem identifier** 104871

C[N+]1(CCC2=CC3=C(C=C2C1C4C5=C(C6=C(C=C5)OCO6)C(=O)O4)OCO3)C.[I-]**SMILES** Source

InChi InChI=1S/C21H20NO6.HI/c1-22(2)6-5-11-7-15-16(26-9-25-15)8-13(11)18(22)19-12-3-4-14-20(27-1

0-24-14)17(12)21(23)28-19; h3-4,7-8,18-19H,5-6,9-10H2,1-2H3; 1H/q+1; /p-1/t18-,19+; /m0./s1-2H3; (H1-q+1)-, /m0./s1-2H3; (H1

InChiKey HKJKCPKPSSVUHY-GRTNUQQKSA-M

MDL number MFCD00078966 **Appearance** Yellow solid

References

Advantages of an antagonist: bicuculline and other GABA antagonists.

Johnston GA (2013) Br J Pharmacol 169(2)

PubMedID 23425285

Differential effects of iontophoretic in vivo application of the GABA(A)-antagonists bicuculline and gabazine in sensory cortex.

Kurt S et al (2006) Hear Res 212(1-2)

PubMedID 16442250

[Bicuculline inhibits airway remodeling in a murine model of chronic asthma].

Zhu T et al (2010) Nan Fang Yi Ke Da Xue Xue Bao 30(4)

PubMedID 20423862