

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

L-AP4

Product overview

Name	L-AP4
Cat No	HB0370
Biological action	Agonist
Purity	>99%
Customer comments	<i>Great product - L-AP4 arrived in time for experiments. Excellent documentation and ease of use. Data looks promising from experiments - happy customer! Verified customer, The University of Newcastle</i>

Description	L-AP4 is a selective group III mGluR agonist which inhibits synaptic transmission.
--------------------	--

Images



Biological Data


Biological description	L-AP4 is a selective group III mGluR agonist (EC_{50} values are 0.9, 252, 0.06-0.6 μM at mGlu₄, mGlu₇, mGlu₈ receptors respectively and >1000 μM at other mGluRs) and also a weak NMDA receptor agonist. L-AP4 is a potent synaptic depressant. L-AP4 reduces glutamate release, inhibits glutamate-mediated EPSPs in the hippocampus, olfactory cortex and spinal cord and inhibits long term potentiation (LTP) in vivo. Additionally, L-AP4 also inhibits GABA release and reversibly reduces GABA-mediated inhibitory post-synaptic potential (IPSP) amplitude. L-AP4 also shows neuroprotective properties.
-------------------------------	---

Solubility & Handling

Storage instructions
Solubility overview
Important

Room temperature
Soluble in water (5mM) and in 0.1M NaOH (100mM)
This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

Chemical Data

Chemical name	L-(+)-2-Amino-4-phosphonobutyric acid
Molecular Weight	183.1
Chemical structure	
Molecular Formula	C ₄ H ₁₀ NO ₅ P
CAS Number	23052-81-5
PubChem identifier	179394
SMILES	C(CP(=O)(O)O)[C@@H](C(=O)O)N
Source	Synthetic
InChi	InChI=1S/C4H10NO5P/c5-3(4(6)7)1-2-11(8,9)10/h3H,1-2,5H2,(H,6,7)(H2,8,9,10)/t3-m/s1
InChiKey	DDOQBQRIEWHWBT-VKHMYYEASA-N
MDL number	MFCD00083244
Appearance	White solid

References

L-AP4, a potent agonist of group III metabotropic glutamate receptor, decreases central action of angiotensin II.

Fedosiewicz-Wasiluk M *et al* (2002) Pol J Pharmacol 54(5)

PubMedID [12593528](#)

The agonist selectivity of a class III metabotropic glutamate receptor, human mGluR4a, is determined by the N-terminal extracellular domain.

Tones MA *et al* (1995) Neuroreport 7(1)

PubMedID [8742431](#)

Changes of mGluR4 and the effects of its specific agonist L-AP4 in a rodent model of diffuse brain injury.

Zhou F *et al* (2003) J Clin Neurosci 10(6)

PubMedID [14592619](#)

Unveiling the functions of presynaptic metabotropic glutamate receptors in the central nervous system.

Schoepp DD (2001) J Pharmacol Exp Ther 299(1)

PubMedID [11561058](#)