

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

CA201019 CellAura fluorescent 5-HT_{2A/2B} receptor antagonist [4F4PP]

Product overview

Name	CA201019 CellAura fluorescent 5-HT _{2A/2B} receptor antagonist [4F4PP]
Cat No	HB7837
Biological description	Fluorescent 5-HT _{2A/2B} receptor antagonist [4F4PP]. Please note that this is an 'in development compound' which has not undergone full biological testing. Please contact us for more information.
Alternative names	CA201019
Biological action	Antagonist
Purity	>97%
Description	Fluorescent 5-HT _{2A} receptor antagonist

Solubility & Handling

Storage instructions	-20 °C (protect from light)
Solubility overview	Soluble in DMSO
Handling	After thawing individual aliquots for use, we recommend briefly sonicating the sample to ensure it is fully dissolved and the solution is homogeneous. We do not recommend using the product after subjecting it to repetitive freeze-thaw cycles.
Shipping conditions	The product, supplied in a dry form, is stable at ambient temperature for periods of up to a few days and does not require shipping on ice/dry ice.
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

Chemical Data

Molecular Weight	987
Source	Synthetic
Formulation	Lyophilized film
Excitation	633 nm
Emission	650 nm
