Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@m2stage.hellobio.com



# **DATASHEET**

Recombinant human CDNF protein

#### **Product overview**

Name Recombinant human CDNF protein

Cat No HB8544

**Biological description** Neurotrophin factor that is expressed in brain, neuronal and some non-neuronal tissues.

CDNF is mainly localized to the endoplasmic reticulum lumen and modulates the unfolded protein

response (UPR) pathway.

It has been shown to promote survival, growth and function of dopamine-specific neurons.

It additionally shows neurorestorative effects in a number of Parkinson's disease animal models. Species of origin

human

**Alternative names** Recombinant Human Cerebral Dopamine Neurotrophic Factor, Cerebral neurotrophic factor, arginine-

rich, mutated in early stage tumors-like 1, Conserved neurotrophic factor, ARMET-like protein 1,

ARMETL1.

**Biological action** Activator **Purity** >96%

Description Neurotrophin implicated in Parkinson's disease research

# **Biological Data**

**Application notes** ED<sub>50</sub> = 15-25μg/ml (determined by its ability to stimulate of proliferation of rat C6 cells) corresponding

to a specific activity of 40-67units/mg.

# **Solubility & Handling**

Solubility overview

To make a stock solution, reconstitute in sterile 18MΩcm water at a concentration > 100µg/ml, which can then be diluted to make a working solution

Handling

- Solutions should be made in sterile deionized water (not less than 100  $\mu g/ml$ ). This solution can then be further diluted with other aqueous solutions.
- Following reconstitution, solutions may be stored at 4°C and are useable for around 2-7 days and for future use store at -18°C.
- For long term storage, a carrier protein (0.1% HSA or BSA) should be added to stock solutions. Solutions should be aliquoted into tightly sealed vials for storage at -20 °C. Freeze-thaw cycles should be prevented.

**Important** 

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

#### **Chemical Data**

**UniProt ID** Q49AH0 E. Coli. Source

White lyophilized powder (sterile filtered & freeze-dried) **Appearance Formulation** Lyophilized from a 0.2µm filtered solution in PBS (pH 7.4)

## References

## **CDNF Protein Therapy in Parkinson's Disease**

Huttunen HJ *et al* (2019) Cell Transplant 28(4) **PubMedID** 30947516

## Characterization of recombinant human brain-derived neurotrophic factor variants

Sunasara KM et al (1999) Arch Biochem Biophys 372(2)

**PubMedID** 10600162

Transport of human recombinant brain-derived neurotrophic factor (BDNF) through the rat blood-brain barrier in vivo using vector-mediated peptide drug delivery

Pardridge WM *et al* (1994) Pharm Res 11(5) **PubMedID**8058646