

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

Recombinant rat CDNF protein

Product overview

Name	Recombinant rat CDNF protein
Cat No	HB8498
Species of origin	rat
Alternative names	Recombinant Rat Cerebral Dopamine Neurotrophic Factor, Cerebral neurotrophic factor, ARMET-like protein 1, Arginine-rich protein mutated in early stage tumors-like 1, Conserved neurotrophic factor, Cdnf, Armetl1.
Purity	>97%
Description	Rat CDNF protein

Biological Data

Application notes	Enhances neurite outgrowth of E16-E18 rat embryonic cortical neurons when immobilized at 5-25 µg/mL on a nitrocellulose-coated microplate.
--------------------------	--

Solubility & Handling

Solubility overview	To make a stock solution, reconstitute in sterile 18MΩcm water at a concentration > 100µg/ml, which can then be diluted to make a working solution
Handling	• Solutions should be made in sterile deionized water (not less than 100 µg/ml). This solution can then be further diluted with other aqueous solutions. • Following reconstitution, solutions may be stored at 4°C and are useable for around 2-7 days and for future use store at -18°C. • For long term storage, a carrier protein (0.1% HSA or BSA) should be added to stock solutions. Solutions should be aliquoted into tightly sealed vials for storage at -20°C. Freeze-thaw cycles should be prevented.
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use

Chemical Data

UniProt ID	P0C5I0
Source	E. Coli.
Appearance	White lyophilized powder (sterile filtered & freeze-dried)
Formulation	Lyophilized from a 0.2µm filtered solution in PBS (pH 7.4)

References

CDNF Protein Therapy in Parkinson's Disease

Huttunen HJ *et al* (2019) Cell Transplant 28(4)
PubMedID [30947516](#)

Characterization of recombinant human brain-derived neurotrophic factor variants

Sunasara KM *et al* (1999) Arch Biochem Biophys 372(2)

PubMedID

10600162

Transport of human recombinant brain-derived neurotrophic factor (BDNF) through the rat blood-brain barrier in vivo using vector-mediated peptide drug delivery

Pardridge WM *et al* (1994) Pharm Res 11(5)

PubMedID

8058646
