Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

TBS-T with Tween 20 (20x)

Product overview

Name TBS-T with Tween 20 (20x)

Cat No HB6971

Biological description 20x Tris Buffered Saline wash buffer for use during immunolabeling techniques such as IHC. When

diluted to 1x contains 130mM NaCl, 20mM Tris, 0.1% Tween 20, 0.002% ProClin-300, pH7.4.

Biological action Buffe

Description 20x Tris Buffered Saline wash buffer to remove non-specific binding and for use during

immunolabeling techniques such as IHC.

Biological Data

Application notes This is a 20x buffer: to use, dilute in dH₂O (50ml of 20x buffer + 950ml dH₂O)

Solubility & Handling

Storage instructions

Important

Room temperature

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

Chemical Data

Appearance Colorless solution