Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

TBS (25x) (pH 7.4)

Product overview

Name TBS (25x) (pH 7.4)

Cat No HB7121

Biological description Tris Buffered Saline (TBS) wash buffer to remove non-specific binding. For use in immunoassays such

as western blot, IHC, ICC, ELISA. Tween 20 may be added to promote more effective washings. The buffer is supplied as a 25X stock solution that should be diluted to 1X with distilled water. When diluted

to 1x contains 130mM NaCl, 20mM Tris, 0.002% ProClin-300, pH7.4.

Biological action

Description Tris Buffered Saline (TBS) wash buffer to remove non-specific binding.

Biological Data

Application notes This is a 25x buffer: to use, dilute in dH_2O (40ml of 20x buffer + 960ml dH_2O)

Solubility & Handling

Storage instructions

Important

Room temperature

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

Chemical Data

Appearance Colorless solution