

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

TBS (25x) (pH 7.4)

Product overview

Name	TBS (25x) (pH 7.4)
Cat No	HB7121
Biological description	Tris Buffered Saline (TBS) wash buffer to remove non-specific binding. For use in immunoassays such as western blot, IHC, ICC, ELISA. Tween 20 may be added to promote more effective washings. The buffer is supplied as a 25X stock solution that should be diluted to 1X with distilled water. When diluted to 1x contains 130mM NaCl, 20mM Tris, 0.002% ProClin-300, pH7.4.
Biological action	Buffer
Description	Tris Buffered Saline (TBS) wash buffer to remove non-specific binding.

Biological Data

Application notes	This is a 25x buffer: to use, dilute in dH ₂ O (40ml of 20x buffer + 960ml dH ₂ O)
--------------------------	--

Solubility & Handling

Storage instructions	Room temperature
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

Chemical Data

Appearance	Colorless solution
-------------------	--------------------
