Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

Product overview

Name Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

Cat No HB7943

Biological description Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue

sections. Optimal for use with primary antibodies that require citrate buffer (pH 6.0) pretreatment. This product contains detergent for emulsification of the paraffin. 1X Dilution: The 100X stock solution

should be diluted 100-fold with distilled water before use.

Biological action ApplicationsBuffer
IHC, IHC-P

Description Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue

sections

Solubility & Handling

Storage instructions

Important

Room temperature

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.