

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

Product overview

Name	Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)
Cat No	HB7943
Biological description	Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue sections. Optimal for use with primary antibodies that require citrate buffer (pH 6.0) pretreatment. This product contains detergent for emulsification of the paraffin. 1X Dilution: The 100X stock solution should be diluted 100-fold with distilled water before use.
Biological action	Buffer
Applications	IHC, IHC-P
Description	Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue sections

Solubility & Handling

Storage instructions	Room temperature
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.
