

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

PBS buffer with Tween 20 (20x)

Product overview

Name	PBS buffer with Tween 20 (20x)
Cat No	HB8088
Biological description	20x wash buffer for use during immunolabeling techniques such as IHC. Tween20 promotes more effective washing to decrease non-specific background staining. When diluted to 1x contains 0.1% Tween-20 in PBS, pH7.4.

Species of origin	Synthetic
Biological action	Buffer
Description	20x wash buffer for use during immunolabeling techniques such as IHC.

Solubility & Handling

Storage instructions	Room temperature
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.
