

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

Citrate Buffer (10x) pH 6.0

Product overview

Name	Citrate Buffer (10x) pH 6.0
Cat No	HB8687
Biological description	Citrate buffer (pH6) for antigen retrieval. Used for antigen retrieval prior to IHC (immunohistochemistry) on formalin-fixed, paraffin-embedded tissue sections mounted on glass slides to increase staining intensity with many primary antibodies.
Biological action	Buffer
Applications	IHC, IHC(IF)
Description	Citrate buffer (pH6) for target retrieval

Solubility & Handling

Storage instructions	Room temperature
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.
