

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com




DATASHEET

Goat Anti-Chicken IgY H&L (HRP) ValidAb™

Product overview

Name	Goat Anti-Chicken IgY H&L (HRP) ValidAb™
Cat No	HB9914
Host	Goat
Clonality	Polyclonal
Target	Chicken IgY H&L
Conjugate	Horseradish peroxidase (HRP)
Description	Goat Anti-Chicken IgG H&L Horseradish peroxidase (HRP) secondary antibody. Part of the ValidAb™ range of highly validated, data-rich antibodies.

Validation data



Product information

Immunogen	Purified chicken IgY
Isotype	IgG
Purification notes	Immunogen affinity chromatography
Concentration	0.6 mg/ml
Formulation	50% Glycerol, 50% PBS pH 7.4 with 0.05% ProClin-300

Tested applications

Applications	ELISA, IHC, WB
Western blot optimal concentration	Dependent on assay. A 1:20,000 dilution is a good starting point.
Open data link	Please follow this link to OSF

Storage & Handling

Storage instructions -20 °C

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use

References

The induction of p53 correlates with defects in the production, but not the levels, of the small ribosomal subunit and stalled large ribosomal subunit biogenesis.

Eastham MJ et al (2023) Nucleic acids research 51

PubMedID [37526268](#)

Western blot: technique, theory, and trouble shooting.

Mahmood T et al (2012) North American journal of medical sciences 4

PubMedID [23050259](#)
