

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@m2stage.hellobio.com



DATASHEET

Goat Anti-Chicken IgY H&L (Janelia Fluor® 646) ValidAb™

Product overview

Name	Goat Anti-Chicken IgY H&L (Janelia Fluor® 646) ValidAb™
Cat No	HB18908
Host	Goat
Clonality	Polyclonal
Target	Chicken IgY H&L
Conjugate	Janelia Fluor® 646
Description	Goat Anti-Chicken IgY H&L Janelia Fluor® 646 secondary antibody. Part of the ValidAb™ range of highly validated, data-rich antibodies.

Product information

Immunogen	Purified chicken IgY
Purification notes	Immunogen affinity chromatography.
Concentration	1 mg/ml
Formulation	20% glycerol in PBS with 0.05% sodium azide and 1% recombinant albumin

Tested applications

Applications	ELISA, FACS and flow cytometry, ICC, live cell imaging, WB, IHC(IF)
IHC(IF) optimal concentration	1:300 to 1:2,000 dilution (0.5 - 3.3 µg/ml). Optimise dependent upon assay. A good starting point is 1:500 (2 µg/ml).
ICC optimal concentration	1:300 to 1:2,000 dilution (0.5 - 3.3 µg/ml). Optimise dependent upon assay. A good starting point is 1:500 (2 µg/ml).

Storage & Handling

Storage instructions	+4 °C
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use
